

**ANTIGENIX AMERCIA Inc.**

**Super-X Plex Flow Cytometry Assay  
Multi-Plex Panel ( 7 Plex)**

Product Information Summary

**Mouse Chemokine      Cat No.MMX370    96 tests**

**DESCRIPTION:** Multi-Plex ( 7 Plex) Flow Cytometry Assay Panel  
- components to measure **Mouse: KC, IP-10, Eotaxin, RANTES,  
MCP-1, MCP-3 and MIG simultaneously-** in Murine Serum,  
Plasma, Cell Culture Supernatant or Cell Lysate samples.

**NOTE: ORDER a Sample Diluent Kit ( to match sample type):**

Mouse/Rat Serum/Plasma (**DKM100**) or  
Cell Culture Supernatant (**DKC100**) or  
Cell Lysate (**DKL100**)

Cell Lysate Prep Buffer (**CLP200**) or Sample Dilution Buffer  
( to dilute high concentration samples) - (**SDL200**)  
also available.

(**STARTER Kit Included**)

**CONTENTS:**

**Antibody Conjugated Beads            5.0 mL (1X) (Pre-Mixed , 7-Plex)**  
**Detection Antibody - Biotin        1.5 mL (2X) (Pre-Mixed , 7-Plex)**  
**Detection Antibody Diluent        1.5 mL (2X)**

**Standard:** One vial of lyophilized standard containing  
**Pre-mixed multiplexed** analytes ( **7-Plex**) ( see vial insert) -  
single standard vial can be used for one or more assays  
indicated.

**STARTER Kit: ( Mouse Panel)**

Reading Buffer            5.0 mL (10X)  
Wash Buffer                20.0 mL (10X)  
Filter plate with lid      ( 96 wells)  
Plate Sealers              (6)  
Streptavidin-PE          3.0 mL (1X)  
PCR 8 tube strip         2 strips

**STORAGE:** Store Refrigerated ( 4-8° C. - in dark)

**CAUTION:** Contains Sodium Azide - ( 0.05%) - Dispose contents  
carefully, flush with water to avoid reaction with lead drain  
lines.

**PROTOCOL:** Refer to Kit protocol provided with assay ordered.

**Assay Specifications:**

	<b>Sensitivities (pg/mL)</b>	<b>Range: (pg/mL)</b>
Mouse KC	< 1.0	2.0 - 2,000
Mouse IP-10	< 3.0	10.0 - 5,000
Mouse Eotaxin	< 3.0	10.0 - 2,000
Mouse RANTES	< 2.0	5.0 - 5,000
Mouse MCP-1	< 1.0	2.0 - 2,000
Mouse MCP-3	< 2.0	5.0 - 5,000
Mouse MIG	< 2.0	5.0 - 5,000

Sample Size: 15.0 uL per test

**X-Reactivity** within panel : **negligable**

Intra Assay CV: <10%      Inter Assay CV <20%

Recovery 70-130% standard dose

**RESEARCH USE ONLY**  
**NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE.**

[WWW.ANTIGENIX.com](http://WWW.ANTIGENIX.com)    1-800-558-1008