

BOSTER BIOLOGICAL TECHNOLOGY Co.,Ltd.

3942 B Valley Ave, Pleasanton, CA, 94566

Phone: 888-466-3604 Fax: 925-215-2184 Email: boster@bosterbio.com Web: www.bosterbio.com

Citrate Buffer

Code AR0024

Size 1 pack (It is sufficient to make 2000 ml of PBS buffer)

Formulation Powder

Content $3g C_6H_8O_7.H_2$, $20g C_6H_5Na_3O_7.2H_2O$

pH Value 6.0

Storage Store at room temperature (RT) for one year

Usage

Reconstitute the powder with distilled water to yield 2000 ml Citrate buffer.

Application

Citrate buffer is used in Heat Induced Epitope Retrieval (HIER) process on formalin-fixed paraffin embedded tissues.