

Anti-Human Interferon Alpha (MMHA-2) Neutralizing Monoclonal Antibody

Catalog Number:	CM2002	Size:	500 µg
Description:	Mouse Anti-Human Interferon Alpha, neutralizing		
Clone:	MMHA-2		
Concentration:	5.0 mg/ml		
Buffer:	phosphate buffered saline (PBS) containing 0.1% bovine serum albumin (BSA)		
Antigen:	human leukocyte interferon		
Isotype:	IgG ₁ κ		
Bioactivity:	neutralizes human interferon alpha; binds to human interferon alpha with high affinity (K _A ≈10 ⁹ M ⁻¹)		

Assay Used to Measure Bioactivity:

One neutralization unit is the amount of antiserum required to neutralize one unit of human interferon alpha (hIFN-α) to a 50% endpoint. Interferon was titrated with the use of the cytopathic effect inhibition assay as described [Rubinstein, S., Familletti, P.C., and Pestka, S. (1981) "Convenient Assay for Interferons," *J. Virol.* **37**, 755-758; Familletti, P.C., Rubinstein, S., and Pestka, S. (1981) "A Convenient and Rapid Cytopathic Effect Inhibition Assay for Interferon," in *Methods in Enzymology*, Vol. 78 (S. Pestka, ed.), Academic Press, New York, 387-394]. In this antiviral assay for interferon about 1 unit/ml of interferon is the quantity necessary to produce a cytopathic effect of 50%. The units are determined with respect to the international reference standard for hIFN-α provided by the National Institutes of Health [see Pestka, S. (1986) "Interferon Standards and General Abbreviations," in *Methods in Enzymology* (S. Pestka, ed.), Academic Press, New York **119**, 14-23].

Usage:	Western blot (1:500); immunoprecipitation; immunohistochemistry; ELISA
Shipping Temperature:	Dry ice
Physical State of Product During Shipping:	Frozen
Special Conditions/Comments:	After receipt, this product should be kept at -20°C or below for retention of full activity. When thawing, the contents of the tube should be apportioned in separate tubes so that freezing and thawing is kept to a minimum.

Usage: For Research Use Only. Not for use in diagnostic procedures.

Country of Origin: USA

* This is a generic data sheet and may be subject to change. Please see the package insert shipped with your product for current data.



Cell Sciences, Inc.
480 Neponset Street
Bldg 12A
Canton, MA 02021

Toll Free: 888-769-1246
Phone: 781-828-0610
Fax: 781-828-0542

E-mail: info@cellsciences.com
Web Site: www.cellsciences.com