



Diacclone

DIACLONE Immunology Products

Anti-Human CD82 Azide Free

Clone: TS82b

CATALOGUE N° :	857.790.000
CLONE :	TS82b
ISOTYPE :	Mouse IgG1
IMMUNISATION :	Jurkat and HEL cell lines
HYBRIDOMA :	Myeloma X63/AG8. x Balb/c spleen cells
SPECIFICITY :	Recognises the CD82 antigen, a 60 kDa protein
QUANTITY :	200 µg / 200 µl
PURIFICATION :	Ion exchange chromatography
FORMAT :	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
STORAGE :	Stable at +2-8°C for 12 months DO NOT FREEZE.
APPLICATION :	Suitable for flow cytometry , immunoprecipitation and Western-Blot (under non-reducing conditions)
REFERENCES :	Bari, R. et al., Am. J. Pathol.,2009; 174(2): 647-660. Charrin S. et al., Eur. J. Immunol., 2003; 33: 2479-2489 Charrin S. et al., Biochem. J., 2003; 373: 409-421 Charrin S. et al., FEBS Lett., 2002; 516: 139-144 Charrin S. et al., J. Biol. Chem., 2001; 276: 14329-14337 Charrin, S. et al., Biochem Biophys Res Commun.,2003; 304(1): 107-12. Gellersen, B. et al., Am. J. Pathol., 2007 ;170(1):126-139 Gellersen, B. et al., Hum. Reprod.,2010; 25(4): 862-873 Gellersen, B., et al., Biol Reprod.,2013; 88(3):80 Gonzalez, M. et al.,Mol. Hum. Reprod.,2011; 17(7):421-433 Poloso, N. J. et al., J Immunol.,2006; 177(8): 5451-8. Ruseva, Z. et al., Exp Cell Res.,2009; 315(10): 1759-71. Xu, C. et al., Faseb J.,2009;fj.08-123414

97.4

58.1

39.8

29

20



1,3 10⁵ HEL cells were lysed (lysis buffer containing 1% Brij97 or triton)
Migration SDS/page 5-15% gel under non-reducing conditions

25 YEARS OF EXPERTISE

With over 25 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

Fast Track Your Research...

DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAplex

DIACLONE SAS

1 Boulevard A. Fleming – BP 1985
25020 Besançon cedex – France
Tel: +33 (0)3 81 41 38 38
Fax: +33 (0)3 81 41 36 36
info@diacclone.com



www.diacclone.com