



Diacclone

## DIACLONE Immunology Products

### Anti-Human CD25 Capture Antibody

Clone: B-G3

- CATALOGUE N° :** 879.500.001  
**CLONE :** B-G3  
**ISOTYPE :** Mouse IgG1  
**IMMUNISATION :** PHA activated T cells  
**HYBRIDOMA :** Myeloma X63/AG.8653 x Balb/c spleen cells  
**SPECIFICITY :** Recognises the Interleukin-2 Receptor alpha chain (IL-2Ra), a 55 kDa protein  
**QUANTITY :** 1.0 mg / 1.0 ml  
**PURIFICATION :** Ion exchange chromatography  
**FORMAT :** Phosphate-buffered saline.  
Sterile-filtered through 0.22 µm  
Carrier and Preservative free  
**STORAGE :** Store at +2-8°C for 12 months  
For longer storage, freeze aliquots at -20°C  
**APPLICATION :** **ELISA Capture Antibody.** This Antibody can be used as Capture in a human CD25 sandwich Immunoassay to detect human CD25 in combination with biotinylated human CD25 Detection Antibody (Cat n° 879.500.002). The **suggested coating concentration** range below should be optimised by each laboratory for each application.

**ELISA:** 0,5-2,5µg/ml

General ELISA protocol is available at  
[www.diacclone.com](http://www.diacclone.com)

#### Related products :

	Contents	Formats	Catalogue Number
CD25 ELISA Set	Capture antibody,	1x96 well	851.650.001
	Biotinylated Detection	5x96 well	851.650.005
	antibody, Standard,	10x96 well	851.650.010
	Streptavidin-HRP, TMB	15x96 well	851.650.015
		20x96 well	851.650.020

Version 02 - 01.11

For research use only

#### 25 YEARS OF EXPERTISE

With over 25 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

#### Fast Track Your Research...

#### DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAPlex

#### DIACLONE SAS

1 Boulevard A. Fleming – BP 1985  
25020 Besançon cedex – France  
Tel: +33 (0)3 81 41 38 38  
Fax: +33 (0)3 81 41 36 36  
[info@diacclone.com](mailto:info@diacclone.com)



[www.diacclone.com](http://www.diacclone.com)