


## EzWay™ DNA 100bp Ladder

|                        |                                                                                                                                                                                                                                                                                                              |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1. Catalog No.         | K0570111                                                                                                                                                                                                                                                                                                     |
| 2. Quantity            | 600 ul                                                                                                                                                                                                                                                                                                       |
| 3. No. of Applications | 120 Assays (5ul per lane)                                                                                                                                                                                                                                                                                    |
| 4. Storage & Stability | -20°C                                                                                                                                                                                                                                                                                                        |
| 5. Description         | EzWay™ DNA 100bp Ladder is ideal for sizing linear double-stranded DNA fragments from 100bp - 15kb, and composed of 12 purified individual DNA fragments. The 500bp and 1kb bands have intensities so they look brighter than other 10 bands. This ladder is pre-mixed with loading dye and is ready to use. |
| 6. Concentration       | Typical bands (500 and 1000bp)<br>Other bands                                                                                                                                                                                                                                                                |
| 7. Contents (bp)       | 12 bands (100, 200, 300, 400, <u>500</u> , 600, 700, 800, 900, <u>1000</u> , 1200, 1500)                                                                                                                                                                                                                     |



2% Agarose Gel

*For research use only; not for use as a diagnostic*

