


EzWay™ DNA 100bp Ladder

1. Catalog No.	K0570115
2. Quantity	600 ul x 5 vials
3. No. of Applications	120 Assays (5ul per lane) x 5 vials
4. Storage & Stability	-20°C
5. Description	EzWay™ DNA 100bp Ladder is ideal for sizing linear double-stranded DNA fragments from 100bp - 15kb, and composed of 12 purified individual DNA fragments. The 500bp and 1kb bands have intensities so they look brighter than other 10 bands. This ladder is pre-mixed with loading dye and is ready to use.
6. Concentration	Typical bands (500 and 1000bp) Other bands
7. Contents (bp)	12 bands (100, 200, 300, 400, <u>500</u> , 600, 700, 800, 900, <u>1000</u> , 1200, 1500)



2% Agarose Gel

For research use only; not for use as a diagnostic



koma@komabiotech.com
Tel. 82-2-2660-5660 82-2-2660-5669