



Anti-Swine fever virus antigen rapid test card

Catalog No. LSY-20010


➤ Principle

This rapid test is an immuno-chromatographic test, using colloidal gold immunoassay method, to detect the antibody against swine fever virus in whole blood, serum and plasma. The whole test only need 20 mins. Easy to operate, fast, accurate, high sensitivity, intuitive result, easy to judge.

➤ Test procedure

- 1) Place the test card on the desktop flatly, mark it;
- 2) Add 2 drops (70~100 ul) serum/ plasma/whole blood sample into the well on the test card;
- 3) Be static at room temperature, then judge the result in 20 minutes.

➤ Results judgement



- 1) Positive: control line and test line are both seen wine red;
- 2) Negative: only control line is seen wine red;
- 3) Invalidation: control line isn't seen wine red.

➤ Notice

- 1) Use the test card as soon as possible once open the aluminum foil, or the test card will be invalidated, store at normal temperature and the shelf life is 18 months.
- 2) If use whole blood as sample, the sample moving speed may be slow, add 1 drop PBS solution 1 minute later after adding sample to help sample moving.

➤ Storage

Specifications: 50 tests/kit

Store at 2~30 °C in sealed dry environment.

Shelf life: 18 months

Shenzhen Lvshiyuan Biotechnology Co., Ltd

D Building, National Biological Industrial Park of Marinelife, No.2 Binhai Road, Dapeng, Shenzhen, 518120 China

Tel. 86-755-28438788

Fax 86-755-28938800

Email: info@lsybt.com

www.lsybt.com
