

Mus musculus NCSTN cDNA Clone



Sino Biological Inc.
Biological Solution Specialist

Catalog Number: MG50934-G

General Information

Gene : nicastrin
Official Symbol : NCSTN
Synonym : Nct, Aph2, AA727311, D1Dau13e, Kiaa0253, mKIAA0253, 9430068N19Rik, Ncstn
Source : *Mus musculus*
cDNA Size: 2127
RefSeq : NM_021607.3

Description

Lot : Please refer to the label on the tube

Sequence Description :

Identical with the Gene Bank Ref. ID sequence except for the point mutations: 62 C/T resulting in the amino acid Ser substitution by Phe, 2033 C/T resulting in the amino acid Thr substitution by Ile, 2038 G/A resulting in the amino acid Val substitution by Ile, 219 T/C, 753 C/T, 1074 T/C and 1609 T/C not causing the amino acid variation.

Vector :

pGEM-T

Shipping carrier :

Each tube contains approximately 10 µg of lyophilized plasmid.

Storage :

The lyophilized plasmid can be stored at ambient temperature for three months.

Quality control :

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list :

M13-47 : 5' GCCAGGGTTTCCAGTCACGAC 3'
RV-M : 5' GAGCGATAACAATTCACACAGG 3'

Other M13 primers can also be used as sequencing primers.

Plasmid Resuspension protocol

1. Centrifuge the tube for 5~10 min at 4,000 rpm.
2. Carefully open the tube and add 100 µl of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than 4000 rpm.
5. Store the plasmid at -20 °C.

The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- *E. coli* transformation
- DNA sequencing

E.coli strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. DH5 α, TOP10, JM109.

Mus musculus NCSTN cDNA Clone



Sino Biological Inc.
Biological Solution Specialist

Catalog Number: MG50934-G

Vector Information

The pGEM-T vector is a high-efficiency TA cloning vector which contains multiple cloning sites as shown below. The pGEM-T vector is 3.0kb in size and contains the ampicin resistance gene for selection. The coding sequence was inserted by TA cloning.

Physical Map of pGEM-T :



pGEM-T Vector



* Please refer to
<http://www.sinobiological.com/Vector-pGEM-T-a-1636.html>
for the vector sequence.

Manufactured By Sino Biological Inc., FOR RESEARCH USE ONLY. NOT FOR USE IN HUMANS.

Fax :+86-10-51029969 • Tel:+86- 400-890-9989 • <http://www.sinobiological.com>