

Product Manual

CytoSelect™ 48- Well Cell Adhesion Assay (ECM Array, Colorimetric Format)

Catalog Number

CBA- 070	48 assays
CBA- 070- 5	5 x 48 assays

FOR RESEARCH USE ONLY
Not for use in diagnostic procedures



CELL BIOLABS, INC.
Creating Solutions for Life Science Research

Introduction

Cell adhesion is a complex process involved in embryogenesis, migration/invasion, tissue remodeling, and wound healing. To perform these processes, cells adhere to extracellular matrix components (via adhesion receptors), forming complexes with components of the cytoskeleton that ultimately affect cell motility, differentiation, proliferation, and survival. The Cell Biolabs CytoSelect™ Cell Adhesion Assay Kit provides a rapid, quantitative method for evaluating cell adhesion. The kit contains sufficient reagents for the evaluation of 48 samples (40 ECM protein-coated wells, 8 BSA-coated wells).

Assay Principle

The CytoSelect™ Cell Adhesion Assay Kit utilizes an ECM protein-coated 48-well plate (see Adhesion Plate Layout below). First, cells are seeded onto the coated substrate, where the adherent cells are captured. Next, unbound cells are washed away, and the adherent cells are fixed/stained. Finally, the stain is extracted and quantified colorimetrically.

Related Products

1. CBA-050: CytoSelect™ 48-Well Cell Adhesion Assay (Fibronectin-Coated, Colorimetric Format)
2. CBA-051: CytoSelect™ 48-Well Cell Adhesion Assay (Fibronectin-Coated, Fluorometric Format)
3. CBA-052: CytoSelect™ 48-Well Cell Adhesion Assay (Collagen I-Coated, Colorimetric Format)
4. CBA-053: CytoSelect™ 48-Well Cell Adhesion Assay (Collagen I-Coated, Fluorometric Format)
5. CBA-056: CytoSelect™ 48-Well Cell Adhesion Assay (Laminin-Coated, Colorimetric Format)
6. CBA-057: CytoSelect™ 48-Well Cell Adhesion Assay (Laminin-Coated, Fluorometric Format)
7. CBA-058: CytoSelect™ 48-Well Cell Adhesion Assay (Fibrinogen-Coated, Colorimetric Format)
8. CBA-059: CytoSelect™ 48-Well Cell Adhesion Assay (Fibrinogen-Coated, Fluorometric Format)
9. CBA-060: CytoSelect™ 48-Well Cell Adhesion Assay (Collagen IV-Coated, Colorimetric Format)
10. CBA-061: CytoSelect™ 48-Well Cell Adhesion Assay (Collagen IV-Coated, Fluorometric Format)
11. CBA-071: CytoSelect™ 48-Well Cell Adhesion Assay (ECM Array, Fluorometric Format)

Kit Components

1. ECM Adhesion Plate (Part No. 107001): One 48-well plate containing 40 ECM protein-coated wells and 8 BSA-coated wells (see layout below). FN, Collagen IV and Fibrinogen are from human, Laminin I is from Mouse and Collagen I is from Bovine.
2. Cell Stain Solution (Part No. 11002): One Bottle – 10.0 mL
3. Extraction Solution (Part No. 11003): One Bottle – 10.0 mL

Adhesion Plate Layout

The following layout indicates the location of wells coated with each ECM protein and those coated with BSA.

	1	2	3	4	5	6	7	8
A	Fibronectin							
B	Collagen I							
C	Collagen IV							
D	Laminin I							
E	Fibrinogen							
F	BSA							

Materials Not Supplied

1. Cell culture medium
2. Serum free medium, such as DMEM containing 0.5% BSA, 2 mM CaCl₂ and 2 mM MgCl₂
3. Cell culture incubator (37°C, 5% CO₂ atmosphere)
4. 1X PBS containing 2 mM CaCl₂ and 2 mM MgCl₂
5. Light microscope
6. 96-well microtiter plate
7. Microtiter plate reader

Storage

Store all kit components at 4°C.

Assay Protocol

1. Under sterile conditions, allow the ECM Adhesion Plate to warm up at room temperature for 10 minutes.
2. Prepare a cell suspension containing 0.1-1.0 x 10⁶ cells/ml in serum free media. Agents that inhibit or stimulate cell adhesion can be added directly to the cell suspension.
3. Add 150 µL of the cell suspension to the inside of each well (BSA-coated wells are provided as a negative control).
4. Incubate for 30-90 min in a cell culture incubator.
5. **Carefully** discard or aspirate the media from each well (**Note: Do not allow wells to dry**). Gently wash each well 4-5 times with 250 µL PBS.

- Aspirate the PBS from each well and add 200 μ L of Cell Stain Solution. Incubate for 10 minutes at room temperature.
- Discard or aspirate the Cell Stain Solution from the wells. Gently wash each well 4-5 times with 500 μ L deionized water.
- Discard the final wash and let the wells air dry.
- Add 200 μ L of Extraction Solution per well, and then incubate 10 minutes on an orbital shaker.
- Transfer 150 μ L from each extracted sample to a 96-well microtiter plate and measure the OD 560nm in a plate reader.

Example of Results

The following figures demonstrate typical results with the CytoSelect™ 48-Well Cell Adhesion Assay Kit. One should use the data below for reference only. This data should not be used to interpret actual results.




Figure 1. ECM-mediated Cell Adhesion. Serum starved cells were allowed to attach to ECM-coated 48-well plate for 1 hr at 100,000 cells/well. Adherent cells were stained (left panel picture) and quantified at OD 560nm after extraction (right panel figure).

References

- Hynes, R. O. (1992) *Cell* **69**:11-25.
- Schwartz, M. A., Schaller, M. D. and Ginsberg, M. H. (1995) *Annu. Rev. Cell Dev. Biol.* **11**:549-599.

Recent Product Citations

- Moore, L.M. et al. (2017). Calcium Sensor, NCS-1, Promotes Tumor Aggressiveness and Predicts Patient Survival. *Mol Cancer Res.* pii: molcanres.0408.2016. doi: 10.1158/1541-7786.MCR-16-0408.

2. Bradbury, P.M. et al. (2017). The focal adhesion targeting (FAT) domain of p130 Crk associated substrate (p130Cas) confers mechanosensing function. *J Cell Sci.* pii: jcs.192930. doi: 10.1242/jcs.192930.
3. Guerra, A.D. et al. (2017). Minocycline enhances the mesenchymal stromal/stem cell pro-healing phenotype in triple antimicrobial-loaded hydrogels. *Acta Biomater.* **51**:184-196. doi: 10.1016/j.actbio.2017.01.021.
4. Ponda, M. and Breslow, J.L. (2016). Serum stimulation of CCR7 chemotaxis due to coagulation factor XIIa-dependent production of high-molecular-weight kininogen domain 5. *PNAS* **113**:E7059-E7068.
5. Grundy, T. J. et al. (2016). Differential response of patient-derived primary glioblastoma cells to environmental stiffness. *Sci Rep.* doi:10.1038/srep23353.
6. Sato, Y. et al. (2016). High mannose-binding *Pseudomonas fluorescens* lectin (PFL) downregulates cell surface integrin/EGFR and induces autophagy in gastric cancer cells. *BMC Cancer.* doi:10.1186/s12885-016-2099-2.
7. Pluciennik, E. et al. (2016). The influence of the WWOX gene on the regulation of biological processes during endometrial carcinogenesis. *Int J Mol Med.* **37**:807-815.
8. Nowakowska, M. et al. (2014). Diverse Effect of WWOX Overexpression in HT29 and SW480 Colon Cancer Cell Lines. *Tumor Biol.* **35**:9291-9301.
9. Jiang, F. et al. (2015). Role of HLA-G1 in trophoblast cell proliferation, adhesion and invasion. *Biochem Biophys Res Commun.* **458**:154-160.
10. Bhat, M. A. et al. (2015). Lead optimization of 2-cyclohexyl-N-[(Z)-(3-methoxyphenyl)/3-hydroxyphenyl] methylidene] hydrazinecarbothioamides for targeting the HER-2 overexpressed breast cancer cell line SKBr-3. *Molecules.* **20**:18246-18263.
11. Pluciennik, E. et al. (2015). The role of WWOX tumor suppressor gene in the regulation of EMT process via regulation of CDH1-ZEB1-VIM expression in endometrial cancer. *Int J Oncol.* **46**:2639-2648.
12. Singla, A. K. et al. (2015). Characterization of a murine model of metastatic human non-small cell lung cancer and effect of CXCR4 inhibition on the growth of metastases. *Oncoscience.* **2**:263-271.
13. Jain, M. et al. (2014). ZNF367 inhibits cancer progression and is targeted by miR-195. *PLoS One.* **9**:e101423.
14. Chanda, D. et al. (2014). Anterior gradient protein-2 is a regulator of cellular adhesion in prostate cancer. *PLoS One.* **9**:e89940.
15. Wang, J. et al. (2014). The iron chelator Dp44mT inhibits hepatocellular carcinoma metastasis via N-Myc downstream-regulated gene 2 (NDRG2)/gp130/STAT3 pathway. *Oncotarget.* **5**:8478.
16. Kucerova, L. et al. (2014) Long-term efficiency of mesenchymal stromal cell-mediated CD-MSC/5FC therapy in human melanoma xenograft model. *Gene Ther.* **21**:874-887.
17. Kośła, K. et al. (2014). WWOX modulates the gene expression profile in the T98G glioblastoma cell line rendering its phenotype less malignant. *Oncol Rep.* **32**:1362-1368.
18. Nowakowska, M, et al. (2014). Diverse effect of WWOX overexpression in HT29 and SW480 colon cancer cell lines. *Tumour Biol.* **35**:9291-9301.
19. Tanoury, Z, et al. (2014). Genes involved in cell adhesion and signaling: a new repertoire of retinoic acid receptor target genes in mouse embryonic fibroblasts. *J. Cell Sci.* **127**:521-533.
20. Xing, C. et al. (2013). Reversing effect of ring finger protein 43 inhibition on malignant phenotypes of human hepatocellular carcinoma. *Mol. Cancer Ther.* **12**:94-103.

21. Mei, J. et al. (2012). Inhibition of IDO1 suppresses cyclooxygenase-2 and matrix metalloproteinase-9 expression and decreases proliferation, adhesion and invasion of endometrial stromal cells. *Mol. Human. Reprod.* 10.1093/molehr/gas021.
22. Kandalam, V. et al. (2011). Lack of tissue inhibitor of metalloproteinases 2 leads to exacerbated left ventricular dysfunction and adverse extracellular matrix remodeling in response to biomechanical stress. *Circulation* **124**:2094-2105.
23. Tsukamoto, H. et al. (2011). Evaluation of anticancer activities of benzo[c]phenanthridine alkaloid sanguinarine in oral squamous cell carcinoma cell line. *Anticancer Res.* **31**:2841-2846.
24. Cervera, A.M. et al. (2008). Cells silenced for SDHB expression display characteristic features of the tumor phenotype. *Cancer Res.* **68**:4058-4067.

Warranty

These products are warranted to perform as described in their labeling and in Cell Biolabs literature when used in accordance with their instructions. THERE ARE NO WARRANTIES THAT EXTEND BEYOND THIS EXPRESSED WARRANTY AND CELL BIOLABS DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR WARRANTY OF FITNESS FOR PARTICULAR PURPOSE. CELL BIOLABS' sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of CELL BIOLABS, to repair or replace the products. In no event shall CELL BIOLABS be liable for any proximate, incidental or consequential damages in connection with the products.

Contact Information

Cell Biolabs, Inc.
7758 Arjons Drive
San Diego, CA 92126
Worldwide: +1 858-271-6500
USA Toll-Free: 1-888-CBL-0505
E-mail: tech@cellbiolabs.com
www.cellbiolabs.com

©2006-2017: Cell Biolabs, Inc. - All rights reserved. No part of these works may be reproduced in any form without permissions in writing.