

RPA150Mu02 50µg

**Recombinant Carcinoembryonic Antigen (CEA)** 

Organism Species: Mus musculus (Mouse)

Instruction manual

FOR RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)



# [PROPERTIES]

Source: Prokaryotic expression

Host: E.coli

Residues: Pro544~Pro771

Tags: N-terminal His Tag

Subcellular Location: n/a

**Purity:** > 97%

Traits: Freeze-dried powder

Buffer formulation: 100mMNaHCO<sub>3</sub>, 500mMNaCl, pH8.3, containing 0.01% SKL, 5%

Trehalose and Proclin300.

Original Concentration: 200µg/mL

**Applications:** Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 6.0

Predicted Molecular Mass: 26.3kDa

Accurate Molecular Mass: 20kDa as determined by SDS-PAGE reducing conditions.

### Phenomenon explanation:

The possible reasons that the actual band size differs from the predicted are as follows:

- 1. Splice variants: Alternative splicing may create different sized proteins from the same gene.
- 2. Relative charge: The composition of amino acids may affects the charge of the protein.
- 3. Post-translational modification: Phosphorylation, glycosylation, methylation etc.
- 4. Post-translation cleavage: Many proteins are synthesized as pro-proteins, and then cleaved to give the active form.
- 5. Polymerization of the target protein: Dimerization, multimerization etc.

# [USAGE]

Reconstitute in 100mM NaHCO3, 500mM NaCl (pH8.3) to a concentration of 0.1-1.0 mg/mL. Do not vortex.



## [STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

# [ SEQUENCE ]

PAYKSHT
SKIVEYHRFT DSATVGAAYR GIEVIFTNGS MVMIDVTEDD AGFYMLEILR
EDFKVEKAYV QLHVNSFVPN SKVSVSTARL SIESVPPSIV GGESVLLLVH
NLPKNLQSLF WYKGVIAEEK SELIQHIIAT SSSLPGPAHS GRETVYSNGS
LLLQRVMQND TGFYTLLTMS TDLKDEIAHV QLQLDTSTCC SLLTSDQLII
VPVPRNIAVG KSVLLLVCNV P

# [ IDENTIFICATION ]




Figure. SDS-PAGE

#### [ IMPORTANT NOTE ]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.