



G-LISA[®] Cdc42 Activation Assay
Biochem Kit[™]
(Absorbance Based)

Cat. # BK127

Manual Contents

Section I: Introduction

Assay Principle -----4

Section II: Kit Contents -----5-6

Section III: Reconstitution and Storage of Components-----7

Section IV: G-LISA[®] Assay

G-LISA[®] Grow Cells and Prepare Lysates -----8-10

G-LISA[®] Assay Preparation-----11

G-LISA[®] Assay Protocol -----12-13

Section V: Technical Support and Purchaser Notification

Troubleshooting -----14

Limited Use Statement -----15

I: Introduction

Assay Principle

The Cdc42 G-LISA[®] kit contains a Cdc42-GTP-binding protein linked to the wells of a 96 well plate. Active, GTP-bound Cdc42 in cell/tissue lysates will bind to the wells while inactive GDP-bound Cdc42 is removed during washing steps. The bound active Cdc42 is detected with a Cdc42 specific antibody. The degree of Cdc42 activation is determined by comparing readings from activated lysates versus non-activated lysates. Inactivation of Cdc42 is generally achieved in tissue culture by a serum starvation step (see Technical Guide). A basic schematic of the steps involved in the G-LISA[®] is shown in Figure 1. Typical G-LISA[®] results are shown in Figure 2.

Figure 1: Simple and Quick Protocol





Figure 2: Typical G-LISA[®] Results



II: Kit Contents

This kit contains enough reagents for 96 assays. When properly stored, kit components are guaranteed stable for a minimum of 6 months. You can assay anywhere from 2 to 96 samples at a time for your own convenience.

Table 1: Kit Contents

Reagents	Cat. # or Part # *	Quantity	Storage
96 well Cdc42-GTP binding plate	Part # GL72	12 strips of 8 wells	Desiccated 4°C
Anti-Cdc42 antibody (mouse monoclonal)	Part # GL14	1 tube, lyophilized	Desiccated 4°C
Secondary antibody - horseradish peroxidase conjugate (HRP)	Part # GL02	2 tubes, lyophilized	Desiccated 4°C
Cdc42 control protein (constitutively active Cdc42)	Part # CDCA	12 tubes, lyophilized	Desiccated 4°C
Cell Lysis Buffer	Part # GL35	1 bottle, lyophilized	Desiccated 4°C
Wash Buffer	Part # PE38	1 tablet	4°C
Antigen Presenting Buffer	Part # GL39	1 bottle, 25 ml	Room temperature
Antibody Dilution Buffer	Part # GL47	1 bottle, lyophilized	Desiccated 4°C
HRP Detection Reagent A	Part # GL43	1 tablet, silver pack	Desiccated 4°C
HRP Detection Reagent B	Part # GL44	1 tablet, gold pack	Desiccated 4°C
HRP Reagent Stop Solution	Part # GL80	2 bottles, 8 ml each	4°C
Precision Red™ Advanced Protein Assay Reagent	Part # GL50 (available as 500 ml size Cat. # ADV02)	1 bottle, 100 ml	Room temperature
Protease Inhibitor Cocktail	Cat. # PIC02	1 tube, lyophilized	Desiccated 4°C

- **Items with part numbers (Part #) are not sold separately and available only in kit format. Items with catalog numbers (Cat. #) are available separately.**

II: Kit Contents Cont'd

The reagents and equipment that you will require but are not supplied:

- Cold 4°C PBS buffer (10 mM phosphate buffer pH 7.4, 140 mM NaCl, 3 mM KCl)
- Concentrated sulfuric acid (need to add 1 ml to HRP Stop Buffer)
- Cell scrapers
- Liquid nitrogen for snap freezing cell lysates.
- Multi-channel or multi-dispensing pipette; 25-200 µl range.
- Multi-channel pipette solution basins (available from VWR Cat. # 21007-970). Used for liquid handling.
- Vortex for mixing cell lysate and other solutions
- Two orbital microplate shakers. Optimal shaker speed is 400 rpm (200 rpm is the minimal speed required). One at room temperature and one at 4°C
- Microplate spectrophotometer (see Technical Guide)

III: Reconstitution and Storage of Components

Many of the kit components are provided in lyophilized form. Prior to beginning the assay, you will need to reconstitute several components as shown in Table 2.

Table 2: Component Storage and Reconstitution

Kit Component	Reconstitution	Storage
96 well Cdc42-GTP binding plate	Keep the plate in the sealed bag with desiccant at all times. Reconstitution is not necessary prior to the start of the assay. If detached, pellets should be tapped to the bottom of the well prior to resuspension.	Store desiccated at 4°C
Anti-Cdc42 antibody	Centrifuge briefly to collect the pellet to the bottom of the tube. Dissolve the powder in 400 µl of PBS.	Store at 4°C
Secondary antibody HRP	Centrifuge briefly to collect the pellet to the bottom of the tube. Dissolve the powder in 80 µl of PBS. <u>Do not use sodium azide</u> in combination with this antibody as it will inactivate the HRP.	Store at 4°C
Cdc42 control protein (12 tubes)	Each tube is good for one experiment. Reconstitution is not necessary until starting the assay (see Table 4).	Store at 4°C
Cell Lysis Buffer	Reconstitute in 100 ml of Milli-Q water. This solution may take 5-10 min to resuspend. Use a 10 ml pipette to thoroughly resuspend the buffer.	Store at 4°C
Wash Buffer	Reconstitute in 1 L of Milli-Q water. This solid will take 45-60 min to resuspend. A magnetic stir bar and stir plate can be used to help resuspension.	Store at room temperature
Antigen Presenting Buffer	No reconstitution necessary.	Store at room temperature
Antibody Dilution Buffer	Reconstitute in 15 ml of Milli-Q water.	Store at 4°C
HRP Detection Reagent A	Resuspend tablet in 10 ml sterile distilled water. Aliquot into 12 x 0.8 ml volumes.	Store at -70°C NOTE: Do not store at -20°C
HRP Detection Reagent B	Resuspend tablet in 10 ml sterile distilled water. Aliquot into 12 x 0.8 ml volumes.	Store at -70°C NOTE: Do not store at -20°C
HRP Stop Solution	Carefully add 1 ml of concentrated sulfuric acid (18 M) to HRP Stop Solution. Check the box on the top of the bottle to indicate acid has been added.	Store at 4°C
Precision Red™ Advanced Protein Assay Reagent	No reconstitution necessary.	Store at room temperature
Protease Inhibitor Cocktail	Reconstitute in 1 ml of dimethyl sulfoxide (DMSO) for 100X stock.	Store at 4°C. The cocktail will freeze at 4°C

IV: How to Grow Cells and Prepare Lysates

First time users should read the G-LISA Technical Guide, especially the description of lysate preparation which is critical for the success of the G-LISA assay. The G-LISA[®] kit uses 50 μ l of lysate (0.25-1 mg/ml lysate protein concentration) per assay. To keep the assay in the linear range, we highly recommend using a lysate concentration of 0.25 - 1 mg/ml for the Cdc42 G-LISA[®].

A) **Growth and Treatment of Cell Lines**

Cells should be plated and grown to desired confluency in appropriate culture conditions (consult literature for particular cell line). Confluent cells can be used experimentally, including for transfection, RNA interference, or serum-starvation, if appropriate. Prior to Cdc42 stimulation, cells should be kept in a “controlled state” via serum starvation so that basal Cdc42 activity is low. Optimal confluency prior to serum starvation and GTPase activation varies by cell line and should be determined empirically. Upon stimulation, Cdc42 proteins are generally activated very rapidly and transiently (30 s to 30 min).

B) **Rapid Processing of Cells to Prepare Lysates**

GTP bound (active) Cdc42 is a labile entity and the bound GTP is susceptible to hydrolysis during and after cell lysis, resulting in Cdc42 inactivation. Rapid processing (<10 min) on ice is essential for accurate and reproducible results. The following guidelines should be followed (See Table 3 for preparing reagents needed for cell lysate preparation).

Washing

1. Retrieve culture dish from incubator, immediately place on ice, aspirate off media, and wash cells with ice-cold PBS to remove serum proteins. NOTE: In some cases Cdc42 can be activated by a cold PBS wash. This should be determined empirically for a specific cell line. In this case room temperature PBS may be a good alternative.
2. Aspirate off all PBS buffer. This is essential so that the lysis buffer is not diluted.

Cell Lysis

To make lysate at a concentration between 0.25 to 1.0 mg/ml, adjust the amount of lysis buffer depending on cell and plate type. Empirically determine the exact lysis volumes for any given cell line. The time period between cell lysis and snap-freezing of lysates is critically important (no more than 10 min on ice). Take the following precautions:

1. Keep solutions and lysates embedded in ice so that the temperature is below 4°C.
2. Lyse cells in an appropriate volume of ice-cold cell lysis buffer.
3. Immediately harvest cell lysates with a cell scraper.
4. Transfer lysates into pre-labeled and pre-chilled 1.5 ml microfuge tubes on ice.
5. Immediately clarify lysates by centrifugation at 10,000 x g , 4°C for 1 min.
6. Save at least 20 μ l of lysate on ice for protein quantification and 50-200 μ l for Cdc42 quantification by western blotting.

IV: How to Grow Cells and Prepare Lysates Cont'd

7. Snap-freeze 120 μ l aliquots of all cell lysates in liquid nitrogen immediately after harvest and clarification. Store at -70°C . Lysates can be stored at -70°C for no longer than 30 days.

C) Measure Lysate Protein Concentration

1. Add 20 μ l of each lysate or lysis buffer into disposable 1 ml cuvettes.
2. Add 1 ml of Precision Red™ Advanced Protein Assay Reagent (Part # GL50) to each cuvette.
3. Incubate for 1 min at room temperature.
4. Blank spectrophotometer with the lysis buffer at 600 nm.
5. Read absorbance of lysates samples.
6. Multiply the absorbance by 5 to obtain the protein concentration in mg/ml (See Technical Guide).
7. Calculate how much ice-cold lysis buffer is needed to equalize the cell extracts to give identical protein concentrations in each sample between 0.25 - 1.0 mg/ml. It is not necessary to equalize protein concentration if the sample variation is less than 10%.

The volume of cold cell lysis buffer to be added to the more concentrated samples can be calculated as follows:

where A is the $\frac{A - B}{B}$ \times (volume of A) = _____ μ l $\frac{\text{higher concentration}}{B}$ is the concentration of the most dilute sample (mg/ml).

Technical Tip

Once an optimal protein concentration for the Cdc42 assay has been determined, it is easier to equalize all cell lysates to the optimal concentration at the time of lysate preparation and before freezing lysate aliquots. This eliminates the need to equalize frozen lysate samples immediately prior to performing the assay and lysates can simply be thawed and used in the G-LISA.

IV: How to Grow Cells and Prepare Lysates Cont'd

Table 3: Reagents Needed for Lysate Preparation

✓	Reagent	Preparation
	Precision Red™ Advanced Protein Assay Reagent	Place on the bench and use at room temperature.
	1.5 ml microfuge tubes, labeled and chilled	Use for aliquoting lysates.
	Ice buckets	Use to pre-chill reagents and scrape cells.
	Protease Inhibitor Cocktail	Resuspend in 1 ml of dimethyl sulfoxide (DMSO) and keep at room temperature.
	Lysis buffer + protease inhibitors, ice-cold	<ol style="list-style-type: none"> Empirically determine volume of Lysis Buffer needed per culture vessel (see Technical Guide). We recommend a final lysate concentration between 0.25-1.0 mg/ml for initial experiments. Determine total volume of Lysis Buffer needed by multiplying the lysis volume per culture vessel (µl) by number of vessels x 1.3. Aliquot this volume of Lysis Buffer into a clearly labeled tube and place in ice. Add 10 µl of protease inhibitor cocktail per ml of aliquoted Lysis Buffer. Mix well and leave on ice. Lysis Buffer needs to be ice cold.
	PBS pH 7.2, ice-cold	Phosphate-buffered saline is not provided in the kit. It should be prepared prior to the assay and placed on ice for at least 30 min to ensure that it is ice cold.
	Cell scrapers	Use to harvest cells.
	Liquid nitrogen	Use to snap-freeze lysate aliquots.

IV: G-LISA Assay Preparation

First time users should read the G-LISA Technical Guide, especially the description of lysate preparation which is critical for the success of the G-LISA assay. The G-LISA[®] kit uses 50 µl of lysate (0.25-1 mg/ml protein concentration) per assay. To keep the assay in the linear range, we highly recommend using a lysate concentration of 0.25 - 1 mg/ml for the Cdc42 G-LISA[®]. The reagents and equipment listed in Table 4 should be prepared prior to performing the assay.

Table 4: Assay Preparation for G-LISA[®]

✓	Reagent	Preparation
	Samples to be assayed	All experimental samples should be prepared prior to G-LISA assay preparation. The following instructions assume that frozen lysates are being used for the assay. Lysates should remain frozen until indicated in G-LISA assay protocol.
	Cdc42-GTP binding 96 well plate	Remove plate from 4°C and keep in its protective bag. Place on your bench at room temperature for 30 min. Do not remove the plate (or strips) from the bag until immediately prior to the experiment.
	Milli-Q water	30 ml placed on ice.
	Protease Inhibitor Cocktail	Resuspend in 1 ml of dimethyl sulfoxide (DMSO) and keep at room temperature.
	Lysis Buffer	Required for blank samples (60 ul per sample) and sample dilutions (if necessary). Add 10 µl of protease inhibitor cocktail per ml of Lysis Buffer. Mix well and leave on ice.
	Cdc42 control protein	Dissolve one tube in 500 µl Lysis Buffer and leave on ice. Use within 15 minutes.
	Anti-Cdc42 antibody	Have primary antibody stock ready on ice. For each 8-well strip, you will need to mix 25 µl antibody with 500 µl Antibody Dilution Buffer (1:20 dilution). This dilution step should be performed just prior to use as detailed in assay protocol.
	Secondary antibody	Have secondary antibody stock ready on ice. For each 8-well strip, you will need to mix 8 µl antibody with 500 µl Antibody Dilution Buffer (1:62.5 dilution). This dilution step should be performed just prior to use as detailed in assay protocol.
	Antibody Dilution Buffer	Place reconstituted buffer on the bench and use at room temperature.
	Wash Buffer	Place on the bench and use at room temperature.
	Antigen Presenting Buffer	Place on the bench and use at room temperature.
	HRP Detection Reagents A and B	The 0.8 ml aliquots of these reagents can remain at -70°C until secondary antibody addition as detailed in the assay protocol.
	HRP Stop Solution	Make sure that the box on top of the bottle is checked, indicating sulfuric acid has been added to the solution. Place the bottle on your bench and allow to warm to room temperature.
	Vortex	Used for mixing reagents, it is helpful to keep one close to the assay area.

IV: G-LISA Assay Protocol

1. Aliquot 120 μ l of Lysis Buffer into a labeled microfuge tube. Place on ice. This is the buffer blank.
2. Mix 30 μ l of CDCA (Cdc42 Control Protein) with 90 μ l of ice-cold Lysis Buffer in a labeled microfuge tube. Place on ice. This is the positive control sample. The positive control is at 2.5ng per well and should be 1-1.3 after background subtraction. Do not re-use.
3. Remove the number of strips required from the Cdc42 plate, place in strip holder, and place on ice. Return remaining strips to storage.
4. Keep the plate on ice and dissolve the powder in the wells with 100 μ l ice-cold water. Detachment of the white powder pellet will not affect assay performance. Tap pellets to the bottom of the wells prior to resuspension.
5. Thaw the snap frozen cell lysates in a ROOM TEMPERATURE water bath. Immediately place on ice after thawing.
6. If not already equalized, add required amount of ice-cold Lysis Buffer to equalize all lysate concentrations. Calculate dilution factors required BEFORE thawing lysates.
7. Completely remove the water from the microplate wells as follows:
Vigorously flick the plate to remove solution from each well, followed by a series of 5-7 vigorous pats onto paper towels. The complete removal of solution from wells between steps of the G-LISA[®] is very important to avoid high background readings. At an absorbance of 490 nm, buffer-only wells should read between 0.10 – 0.4 and positive control wells should read between 0.4-0.6 (after subtraction of blank).
8. Return plate back to ice. Immediately add 50 μ l of equalized cell lysate to replicate wells.
9. Pipette 50 μ l of buffer blank control into duplicate wells.
10. Pipette 50 μ l of Cdc42 positive control into duplicate wells.
11. Immediately place the plate on a cold orbital microplate shaker (400 rpm recommended, 200 rpm minimum) at 4°C for exactly 15 min.
12. During the incubation, dilute the anti-Cdc42 primary antibody to 1/20 in Antibody Dilution Buffer (add 25 μ l of antibody to every 500 μ l Antibody Dilution Buffer). Note: The final volume of 500 μ l is adequate for one strip (8 wells).
13. After 15 min, remove the solution from the wells and wash twice with 200 μ l Wash Buffer at room temperature using a multi-channel pipettor. Do not leave plate unattended at this time. Vigorously remove the Wash Buffer after each wash as described in step 7.
14. Place plate on the bench.
15. **Immediately pipette 200 μ l of room temperature Antigen Presenting Buffer into each well using a multi-channel pipettor and incubate at room temperature for exactly 2 min.**

IV: Assay Protocol Cont'd

16. Vigorously flick out the Antigen Presenting Buffer as described in step 7.
16. Immediately wash the wells three times with 200 μ l of room temperature Wash Buffer, removing Wash Buffer as described in step 7.
17. Add 50 μ l of diluted anti-Cdc42 primary antibody to each well and leave the plate on the orbital microplate shaker (200-400 rpm) at room temperature for 30 min.
18. During primary antibody incubation, dilute secondary HRP labeled antibody to 1/62.5 in Antibody Dilution Buffer (add 8 μ l of antibody to every 500 μ l Antibody Dilution Buffer). Note: The final volume of 500 μ l is adequate for one strip (8 wells).
19. Vigorously flick out the anti-Cdc42 primary antibody as described in step 7.
20. Immediately wash the wells three times with 200 μ l of room temperature Wash Buffer, removing Wash Buffer as described in step 7.
21. Add 50 μ l of diluted Secondary antibody to each well and leave the plate on a microplate shaker (200–400 rpm) at room temperature for 30 min.
22. During secondary antibody incubation, thaw an aliquot of HRP detection reagents A and B in a room temperature water bath and remove as soon as they are thawed. Do not mix.
23. Immediately prior to the end of the secondary antibody incubation, mix HRP detection reagents A and B in equal volumes (70 μ l of A/B mixture per well is needed). Protect mixture from light. Discard unused solution.
24. Vigorously flick out the secondary antibody as described in step 7.
25. Wash the wells three times with 200 μ l of room temperature Wash Buffer as described in step 7.
26. Pipette 70 μ l of the mixed HRP detection reagent into each well and incubate at 37°C for 15 min.
27. Add 140 μ l of HRP Stop Buffer to each well.
28. Check that the wells are free of bubbles; if not, remove before continuing.
29. Immediately read the signal by measuring absorbance at 490 nm using a microplate spectrophotometer. Designate Lysis Buffer only wells as the assay Blank.

V: Troubleshooting

Observation	Possible cause	Remedy
Weak or no signal in all wells.	<ol style="list-style-type: none"> 1. Slow processing of samples or processing > 4°C. 2. Wells dried out during experiment. 3. Plate became damp during storage. Well contents will appear sticky and opaque. 4. A step or component of the assay was omitted. 5. Insufficient HRP reaction time. 	<ol style="list-style-type: none"> 1. Process samples quickly on ice. Snap freeze aliquots. 2. Do not remove the solution in the wells unless the next solution is ready. 3. Store the plate in the desiccant bag with the bag securely sealed. Keep the cover on the plate. If wells appear sticky and opaque, the plate can no longer be used. 4. Confirm with checklist that all reagents were added. 5. When Cdc42 signal is very low, allow longer HRP development times for a stronger signal. HRP Stop Solution should be added prior to reading at 490 nm.
High signal in all wells.	<ol style="list-style-type: none"> 1. Antibody concentration is too high. 2. Washes were insufficient. 3. HRP reagents stored incorrectly 	<ol style="list-style-type: none"> 1. Follow the recommended dilution of antibodies in the manual. If still too high, an antibody titration is necessary to optimize your results. 2. Follow the instructions for washing in the manual. 3. Aliquot HRP solutions and store at -70° C. Do not reuse aliquots once thawed.
Background readings are high (>0.40).	<ol style="list-style-type: none"> 1. Inefficient removal of solutions from G-LISA® wells. 2. HRP reagents stored incorrectly 	<ol style="list-style-type: none"> 1. Background should read between 0.10 – 0.40. Vigorous flicking and patting of the inverted plate is required to <u>completely</u> remove solutions from the wells after each step is complete. See G-LISA instructional video for details. 2. Aliquot HRP solutions and store at -70° C. Do not reuse aliquots once thawed
Induced sample does not give significant signal increase.	<ol style="list-style-type: none"> 1. Poor inducer activity. 2. Technique not rapid or cold enough. 4. The endogenous GTP-Cdc42 level is too high. 5. Tissue culture cells not correctly serum starved. 6. Temperature of lysis and incubation is not 4°C. 	<ol style="list-style-type: none"> 1. Purchase a fresh vial of inducer. 2. Confirm instructions were followed using the Experiment Record Sheet (see Technical Guide). 4. Titrate down the amount of extract to be added. 5. See Technical Guide and referenes therein for guidance on serum starvation. It is a good idea to stain cells with phalloidin to qualitatively determine success of serum starvation and induction. 6. Lyse cells on ice, keep cell lysis buffer and distilled water on ice.
Significant variation between replicate samples.	<ol style="list-style-type: none"> 1. Incorrect volume of solutions for each step added in the wells. 2. Inaccurate pipetting. 	<ol style="list-style-type: none"> 1. Follow the instruction for recommended volume in the manual. 2. A multi-channel pipettor is recommended.
Positive control not working.	<ol style="list-style-type: none"> 1. Positive control protein was stored incorrectly after reconstitution. 	<ol style="list-style-type: none"> 1. Each tube of positive control should be used within 15 minutes. Discard any unused control protein as it is not stable.

V: Limited Use Statement

Limited Use Statement

The G-LISA[®] kits are based on patented technology developed at Cytoskeleton Inc. (Patent# 7,763,418 B2). The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of product and components of product in research conducted by the buyer. The buyer cannot sell or otherwise transfer this product or any component thereof to a third party or otherwise use this product or its components for commercial purposes. Commercial purposes include, but are not limited to: use of the product or its components in manufacturing; use of the product or its components to provide a service; resale of the product or its components.

The terms of this Limited Use Statement apply to all buyers including academic and for-profit entities. If the purchaser is not willing to accept the conditions of this Limited Use Statement, Cytoskeleton Inc. is willing to accept return of the unused product with a full refund.

