

Human Recombinant AMY1 Receptor Stable Cell Line

Cat. No. M00557

Version 03232016

I	INTRODUCTION	1
II	BACKGROUND.....	1
III	REPRESENTATIVE DATA.....	2
IV	THAWING AND SUBCULTURING.....	2
V	REFERENCES.....	3
	Limited Use License Agreement.....	4

I. INTRODUCTION

Catalog Number: M00557

Cell Line Name: CHO-K1/G α 15/AMY1

Gene Synonyms: RAMP1+CALCR

Expressed Gene: Genbank Accession Number NM_005855+NM_001742; no expressed tags

Host Cell: CHO-K1/G α 15

Quantity: Two vials of frozen cells (3 \times 10⁶ per vial)

Stability: 16 passages

Application: Functional assay for AMY1 receptor

Freeze Medium: 45% culture medium, 45% FBS, 10% DMSO

Complete Growth Medium: Ham's F12, 10% FBS

Culture Medium: Ham's F12, 10% FBS, 200 μ g/ml Zeocin, 100 μ g/ml Hygromycin B

Mycoplasma Status: Negative

Storage: Liquid nitrogen immediately upon delivery

II. BACKGROUND

Receptor activity-modifying proteins (RAMPs) are a class of protein which interact with and modulate the activities of several Class B G Protein-Coupled Receptors including the receptors for secretin, calcitonin (CT), glucagon, and vasoactive intestinal peptide (VIP). There are three distinct types of RAMPs, designated RAMP1, RAMP2, and RAMP3, each encoded by a separate gene.

Currently the function of RAMPs is divided into 2 class activities. Association of RAMPs with either the CT or CALCRL proteins forms 6 different receptors from the calcitonin receptor family. When associated with the Calcitonin receptor (CTR) or Calcitonin receptor-like (CALCRL), RAMPs can change the selectivity of the receptor for a specific hormone. In the cases of the other receptors mentioned however, there is no evidence that they can do this, but instead function to regulate trafficking of receptors from the ER / golgi to the membrane.

GenScript's cloned human AMY1-expressing cell line is generated in the CHO-K1/G α 15 host.

§: GenScript employs a PCR-based method to test the mycoplasma. The test covers 11 of the most common strains of mycoplasma, (covering approximately 95% of *M. fermentans*, *M. hyorhinis*, *M. arginini*, *M. orale*, *M. salivarium*, *M. hominis*, *M. pulmonis*, *M. arthritidis*, *M. neurolyticum*, *M. hyopneumoniae* and *M. capricolum*) and one species *Ureaplasma (U. urealyticum)*, with sufficient sensitivity and specificity.

860 Centennial Ave., Piscataway, NJ 08854, USA

Toll-Free: 1-877-436-7274 Tel: 1-732-885-9188 Fax: 1-732-210-0262 Email: product@genscript.com Web: www.genscript.com

III. REPRESENTATIVE DATA

Concentration-dependent stimulation of intracellular calcium mobilization by α -CGRP in CHO-K1/ $G_{\alpha 15}$ /AMY1 and CHO-K1/ $G_{\alpha 15}$ cells




Figure 1. α -CGRP-induced concentration-dependent stimulation of intracellular calcium mobilization in CHO-K1/ $G_{\alpha 15}$ /AMY1 and CHO-K1/ $G_{\alpha 15}$ cells. The cells were loaded with Calcium-4 prior to stimulation with an AMY1 receptor agonist, α -CGRP. The intracellular calcium change was measured by FlexStation. The stimulation% were plotted against the log of the cumulative doses (5-fold dilution) of α -CGRP (Mean \pm SD, n = 4). The EC₅₀ of α -CGRP on AMY1 expressing in CHO-K1/ $G_{\alpha 15}$ cells was 58.8 nM.

Notes:

- EC₅₀ value is calculated with four parameter logistic equation:

$$Y = \text{Bottom} + (\text{Top} - \text{Bottom}) / (1 + 10^{((\text{LogEC}_{50} - X) * \text{HillSlope}))}$$

X is the logarithm of concentration. Y is the response
Y is RFU and starts at Bottom and goes to Top with a sigmoid shape.
- Signal to background Ratio (S/B) = Top/Bottom

IV. THAWING AND SUBCULTURING

Thawing Protocol

- Remove the vial from liquid nitrogen tank and thaw cells quickly in a 37°C water-bath.
- Just before the cells are completely thawed, decontaminate the outside of the vial with 70% ethanol and transfer the cells to a 15 ml centrifuge tube containing 9 ml of complete growth medium.
- Pellet cells by centrifugation at 200 x g force for 5 min, and remove the medium.
- Resuspend the cells in complete growth medium.
- Transfer the cell suspension to a 10 cm dish with 10 ml of complete growth medium.
- Grow the cells in incubator with 37°C, 5 %CO₂.
- In the following day, replace the cells with fresh medium contains antibiotic.

860 Centennial Ave., Piscataway, NJ 08854, USA

Toll-Free: 1-877-436-7274 Tel: 1-732-885-9188 Fax: 1-732-210-0262 Email: product@genscript.com Web: www.genscript.com

Sub-culturing Protocol

1. Remove the culture medium from cells.
2. Wash cells with PBS (pH=7.4) to remove all traces of serum that contains trypsin inhibitor.
3. Add 2.0 ml of 0.05% (w/v) Trypsin- EDTA (GIBCO, Cat No. 25300) solution into 10 cm dish and observe the cells under an inverted microscope until cell layer is dispersed (usually within 3 to 5 minutes).
Note: To avoid cells clumping, do not agitate the cells by hitting or shaking the dish while waiting for the cells detach. If cells are difficult to detach, please place the dish in 37°C incubator for ~2 min.
4. Add 6.0 to 8.0 ml of complete growth medium into dish and aspirate cells by gently pipetting.
5. Centrifuge the cells at 200 x g force for 5min, and remove the medium.
6. Re-suspend the cells in culture medium and add the cells suspension to new culture dish.
7. Grow the cells in incubator with 37°C, 5 %CO₂.

Subcultivation Ratio: 1:3 to 1:8 weekly.

Medium Renewal: Every 2 to 3 days

V. REFERENCES

1. Sexton PM, Morfis M, Tilakaratne N, et.al(2006). Complexing receptor pharmacology: modulation of family B G protein-coupled receptor function by RAMPs. *Ann N Y Acad Sci.* 1070: 90–104.
2. Young A (2005). Amylin: Physiology and Pharmacology; Chapter 3: Receptor pharmacology. *Advances in Pharmacology* 52: 47–65.
3. *Calcitonin Receptors: Introduction. IUPHAR Database of Receptors and Ion Channels. International Union of Basic and Clinical Pharmacology.
4. McLatchie LM, Fraser NJ, Main MJ, et. al(1998). RAMPs regulate the transport and ligand specificity of the calcitonin-receptor-like receptor. *Nature* 393 (6683): 333–9.
5. Foord SM, Marshall FH (1999). RAMPs: accessory proteins for seven transmembrane domain receptors. *Trends Pharmacol. Sci.* 20 (5): 184–7.

GenScript USA Inc,
860 Centennial Ave.
Piscataway, NJ 08854
Toll-Free: 1-877-436-7274
Tel: 1-732-885-9188, Fax: 1-732-210-0262
Email: product@genscript.com
Web: <http://www.genscript.com>

For Research Use Only.

860 Centennial Ave., Piscataway, NJ 08854, USA

Toll-Free: 1-877-436-7274 Tel: 1-732-885-9188 Fax: 1-732-210-0262 Email: product@genscript.com Web: www.genscript.com

Limited Use License Agreement

This is a legal agreement between you (Licensee) and GenScript USA Inc. governing use of GenScript's stable cell line products and protocols provided to licensee. By purchasing and using the stable cell line, the buyer agrees to comply with the following terms and conditions of this label license and recognizes and agrees to such restrictions:

- 1) The products are not transferable and will be used at the site where they were purchased. Transfer to another site owned by buyer will be permitted only upon written request by buyer followed by subsequent written approval by GenScript.
- 2) The purchaser cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party.

The products sold by GenScript are for laboratory and animal research purposes only. The products are not to be used on humans, for consumption, or for any unlawful uses.

GenScript USA Inc. will not assert against the buyer a claim of infringement of patents owned or controlled by GenScript USA Inc. and claiming this product based upon the manufacture, use or sale of a clinical diagnostic, therapeutic and vaccine, or prophylactic product developed in research by the buyer in which this product or its components has been employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on the use of this product for other purposes, contact Marketing Department, GenScript USA Inc., 120 Centennial Avenue, Piscataway, New Jersey 08840, U.S.A. Phone: 1-732-885-9188. Fax: 1-732-210-0262. Email: marketing@genscript.com.