


Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK  
Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: [info@abbexa.com](mailto:info@abbexa.com)

### 1Kb Plus DNA Marker

Catalogue No.: abx098051



1 Kb Plus DNA Marker is composed of twelve linear double-stranded DNA bands. This ready-to-use product is premixed with 1xDNA loading buffer. Based on experiment requirement, 5 µl marker can be directly loaded for gel electrophoresis. It provides simple operation and produces clear image.

**Target:** 1Kb Plus DNA Marker

**Storage:** Store at -20°C. Stable for 12 months from date of receipt.

**Buffer:** 10 mM Tris-HCl (pH 8.4), 10 mM EDTA, 0.02% bromophenol blue, 5% glycerol.

**Note:** This product is for research use only.