


Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK
Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

100bp Plus DNA Marker (12 linear double-stranded DNA bands)

Catalogue No.: abx098053

100bp Plus DNA Marker abx098053



100bp Plus DNA Marker is comprised of twelve linear double-stranded DNA bands. This ready-to-use product is premixed with 1x DNA loading buffer. Based on experiment requirement, 5 µl marker can be directly loaded for gel electrophoresis. It provides simple operation and produces clear image.

Target: 100bp Plus DNA Marker

Storage: Store at -20°C. Stable for 12 months from date of receipt.

Buffer: 10 mM Tris-HCl (pH 8.4), 10 mM EDTA, 0.02% bromophenol blue, 5% glycerol.

Note: This product is for research use only.