


Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK
Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

25-90 kDa Protein Marker

Catalogue No.: abx098119



Protein Marker is composed of eight proteins ranging from 25 kDa to 90 kDa. Protein bands of 25 kDa, 45 kDa and 85 kDa are prestained allowing easy identification and monitoring of your proteins during electrophoresis and membrane transfer. Other five proteins contain several IgG binding sites, allowing marker visualization using the same reagents and protocol for your target proteins. These no-dye attached proteins result in more accurate molecular weight estimation. Prestained bands for monitoring electrophoresis and membrane transfer. Visualizing molecular weight marker on western blots with the commonly used immunodetection methods. Ready-to-use format, no heating and reductant are required.

Target: 25-90 kDa Protein Marker

Storage: Store at -20°C. Stable for 12 months from date of receipt.

Buffer: 125 mM Tris-HCl (pH 6.8), 5 mM EDTA, 10 mM DTT, 17.4% glycerol, 3% SDS, 0.01% phenol red.

Note: This product is for research use only.