


## **DATASHEET**

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

## 30-150 kDa Protein Marker

Catalogue No.:abx098121



30-150 kDa Protein Marker is composed of seven proteins ranging from 30 kDa to 150 kDa. Protein bands of 30 kDa and 120 kDa are prestained allowing easy identification and monitoring of your proteins during electrophoresis and membrane transfer. Other five proteins contain several IgG binding sites, allowing marker visualization using the same reagents and protocol for your target proteins. These no-dye attached proteins result in more accurate molecular weight estimation. Prestained bands for monitoring electrophoresis and membrane transfer. Visualizing molecular weight marker on western blots with the commonly used immuno-detection methods. Ready-to-use format, no heating and reductant are required.

Target: 30-150 kDa Protein Marker

**Storage:** Store at -20°C. Stable for 12 months from date of receipt.

Buffer: 125 mM Tris-HCl (pH 6.8), 5 mM EDTA, 10 mM DTT, 17.4% glycerol, 3% SDS, 0.01% phenol red.

**Note:** This product is for research use only.