

SARS-COV-2 Spike S Trimer Protein

Cat. No. COV-VM4SS



Description

Source	Recombinant SARS-COV-2 Spike S is expressed from HEK293 with His tag and Avi tag at the C-Terminus. It contains Val16-Glu1188.
Accession	QHD43416.1
Molecular Weight	The protein has a predicted MW of 136.6 kDa. Due to glycosylation, the protein migrates to 175-230 kDa based on Tris-Bis PAGE result.
Endotoxin	Less than 1EU per µg by the LAL method.
Purity	> 95% as determined by Tris-Bis PAGE

Formulation and Storage


Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -20 to -80°C for 3-6 months in unopened state after reconstitution. 2-8°C for 2-7 days after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

The spike protein (S) of coronavirus (CoV) attaches the virus to its cellular receptor, angiotensin-converting enzyme 2 (ACE2). A defined receptor-binding domain (RBD) on S mediates this interaction. The S protein plays key parts in the induction of neutralizing-antibody and T-cell responses, as well as protective immunity.

Assay Data

Tris-Bis PAGE




SARS-COV-2 Spike S Trimer on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.

ELISA Data

Assay Data

SARS-COV-2 Spike S Trimer, His Tag ELISA

0.1µg SARS-COV-2 Spike S Trimer, His Tag Per Well



Immobilized SARS-COV-2 Spike S Trimer, His Tag at 1µg/ml (100µl/Well). Dose response curve for Human ACE2, hFc Tag with the EC50 of 42.2ng/ml determined by ELISA.