

SARS-CoV Spike RBD Antibody, Rabbit PAb, Antigen Affinity Purified

Has cross-reactivity in ELISA with 2019-nCoV Spike RBD (Cat# 40592-V05H), 2019-nCoV Spike S1 (Cat# 40591-V08H).



Sino Biological
Biological Solution Specialist

Catalog Number: 40150-T52

GENERAL INFORMATION

Immunogen:	Recombinant SARS-CoV Spike RBD Protein (Catalog#40150-V08B2)
Preparation	Produced in rabbits immunized with purified, recombinant SARS-CoV Spike RBD (Catalog#40150-V08B2; AAX16192.1; Arg306-Phe527). SARS-CoV Spike RBD specific IgG was purified by SARS-CoV Spike RBD affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	SARS-CoV Spike RBD Protein Has cross-reactivity in ELISA with 2019-nCoV Spike RBD (Cat# 40592-V05H), 2019-nCoV Spike S1 (Cat# 40591-V08H).
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

APPLICATIONS

Applications:	WB,ELISA,FCM,IP
	(Antibody's applications have not been validated with corresponding viruses. Optimal concentrations/dilutions should be determined by the end user.)

RECOMMENDED CONCENTRATION

Western Blot	WB: 1:1000-1:5000
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human SARS Coronavirus Spike RBD .

Please Note: Optimal concentrations/dilutions should be determined by the end user.


SARS-CoV Spike RBD Antibody, Rabbit PAb, Antigen Affinity Purified

Cross-reactivity in ELISA with 2019-nCoV Spike RBD (Cat# 40592-V05H),
2019-nCoV Spike S1 (Cat# 40591-V08H)



Sino Biological
Biological Solution Specialist

Catalog Number: 40150-T52



Anti-Human SARS Coronavirus Spike rabbit polyclonal antibody at 1:1000 dilution.
Sample: Human SARS Coronavirus Spike Recombinant Protein
Lane A: 50ng
Lane B: 10ng

Secondary
Goat Anti- Rabbit IgG H&L (Dylight 800) at 1/10000 dilution.

Developed using the Odyssey technique.
Performed under reducing conditions.